

For Research Use Only

TRMT12 Recombinant antibody

Catalog Number: 84167-3-RR



Basic Information

Catalog Number:

84167-3-RR

Size:

100ul , Concentration: 1000 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9258

GenBank Accession Number:

BC011713

GeneID (NCBI):

55039

UNIPROT ID:

Q53H54

Full Name:

tRNA methyltransferase 12 homolog (S. cerevisiae)

Calculated MW:

448 aa, 50 kDa

Observed MW:

45-50 kDa

Purification Method:

Protein A purification

CloneNo.:

241483H1

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells, K-562 cells, MCF-7 cells, U2OS cells, U-521 cells, MDA-MB-231 cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84167-3-RR (TRMT12 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.