

For Research Use Only

# PRKAG2 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 84169-5-PBS



## Basic Information

<b>Catalog Number:</b> 84169-5-PBS	<b>GenBank Accession Number:</b> BC020540	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ug , Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 51422	<b>CloneNo.:</b> 260578C10
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q9UGJ0	
<b>Isotype:</b> IgG	<b>Full Name:</b> protein kinase, AMP-activated, gamma 2 non-catalytic subunit	
<b>Immunogen Catalog Number:</b> AG3214	<b>Calculated MW:</b> 569 aa, 63 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Product Information

84169-5-PBS targets PRKAG2 as part of a matched antibody pair:

MP01102-2: 84169-6-PBS capture and 84169-5-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage

**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



The mean PRKAG2 concentration was determined to be 26.9 ng/mL in HeLa cell extract based on a 1.0 mg/mL extract load, 1.3 ng/mL in HCT116 cell extract based on a 1.0 mg/mL extract load.



Cytometric bead array standard curve of MP01102-2, PRKAG2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84169-6-PBS. Detection antibody: 84169-5-PBS. Standard: Ag3214. Range: 0.625-80 ng/mL.