

For Research Use Only

Ki-67 Recombinant antibody

Catalog Number: 84192-5-RR



Basic Information

Catalog Number: 84192-5-RR	GenBank Accession Number: NM_002417	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 4288	CloneNo.: 241499C6
Source: Rabbit	UNIPROT ID: P46013	Recommended Dilutions: IHC: 1:250-1:1000
Isotype: IgG	Full Name: antigen identified by monoclonal antibody Ki-67	
	Calculated MW: 359 kDa	

Applications

Tested Applications: IHC, ELISA	Positive Controls: IHC : human tonsillitis tissue,
Species Specificity: human	
Note-IHC: suggested antigen retrieval with <i>TE buffer pH 9.0</i>; (*) Alternatively, antigen retrieval may be performed with <i>citrate buffer pH 6.0</i>	

Storage

Storage: Store at -20°C. Stable for one year after shipment.
Storage Buffer: PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20° C storage

*** 20ul sizes contain 0.1%BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 84192-5-RR (KI67 antibody) at dilution of 1:500 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).