

For Research Use Only

ANKRD6 Recombinant antibody, PBS Only

Catalog Number: 84210-4-PBS



Basic Information

Catalog Number:

84210-4-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19328

GenBank Accession Number:

BC042173

GeneID (NCBI):

22881

UNIPROT ID:

Q9Y2G4

Full Name:

ankyrin repeat domain 6

Calculated MW:

727 aa, 80 kDa

Purification Method:

Protein A purification

CloneNo.:

241494D11

Applications

Tested Applications:

IF/ICC, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed A549 cells using ANKRD6 antibody (84210-4-RR, Clone: 241494D11) at dilution of 1:250 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 84210-4-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 84210-4-RR against Human ANKRD6 were performed. The affinity constant is 0.722 nM.