

For Research Use Only

ELL2 Recombinant monoclonal antibody, PBS Only

Catalog Number: 84226-4-PBS



Basic Information

Catalog Number: 84226-4-PBS	GenBank Accession Number: BC028412	Purification Method: Protein A purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 22936	CloneNo.: 260362A10
Source: Rabbit	UNIPROT ID: O00472	
Isotype: IgG	Full Name: elongation factor, RNA polymerase II, 2	
Immunogen Catalog Number: AG3428	Calculated MW: 640 aa, 72 kDa	
	Observed MW: 80-90 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

ELL2 (Elongation Factor for RNA Polymerase II 2) is a key transcriptional elongation factor that regulates RNA polymerase II (Pol II) during mRNA synthesis (PMID: 34265249). It belongs to the ELL gene family, which includes ELL1, ELL2, and ELL3, and plays a crucial role in controlling the elongation phase of transcription (PMID: 9108030). ELL2 is an androgen-responsive gene that is expressed by prostate epithelial cells and is frequently down-regulated in prostate cancer (PMID: 29179998). The predicted molecular weight of the ELL2 protein is 72 kDa, but due to post-translational modifications, its molecular weight is approximately 90 kDa. (PMID: 39582094, PMID: 34948391)

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84226-4-RR (ELL2 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 84226-4-PBS in a different storage buffer formulation.



WB result of ELL2 antibody (84226-4-RR; 1:400; incubated at room temperature for 1.5 hours) with sh-Control and sh-ELL2 transfected HeLa cells. This data was developed using the same antibody clone with 84226-4-PBS in a different storage buffer formulation.