

For Research Use Only

# Beta TrCP Recombinant antibody, PBS Only

Catalog Number: 84297-1-PBS



## Basic Information

Catalog Number:

84297-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG28739

GenBank Accession Number:

BC027994

GeneID (NCBI):

8945

UNIPROT ID:

Q9Y297

Full Name:

beta-transducin repeat containing

Calculated MW:

605 aa, 69 kDa

Purification Method:

Protein A purification

CloneNo.:

241287B8

## Applications

Tested Applications:

IF/ICC, Indirect ELISA

Species Specificity:

human

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed MCF-7 cells using Beta TrCP antibody (84297-1-RR, Clone: 241287B8 ) at dilution of 1:250 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 84297-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 84297-1-RR against Human Beta TrCP were performed. The affinity constant is 2.81 nM.