

For Research Use Only

SKOR2 Recombinant antibody

Catalog Number: 84303-1-RR



Basic Information

Catalog Number:

84303-1-RR

Size:

100ul , Concentration: 1000 µg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

XM_001715219

GeneID (NCBI):

652991

UNIPROT ID:

Q2VWA4

Full Name:

functional smad suppressing element 18

Calculated MW:

104 kDa

Observed MW:

65 kDa

Purification Method:

Protein A purification

CloneNo.:

241181A6

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : SH-SY5Y cells, mouse brain tissue, mouse cerebellum tissue, rat brain tissue, rat cerebellum tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84303-1-RR (SKOR2 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.