

For Research Use Only

POLR3C Recombinant antibody, PBS Only

Catalog Number: 84367-7-PBS



Basic Information

Catalog Number:

84367-7-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27311

GenBank Accession Number:

BC002586

GeneID (NCBI):

10623

UNIPROT ID:

Q9BUI4

Full Name:

polymerase (RNA) III (DNA directed)
polypeptide C (62kD)

Calculated MW:

61 kDa

Purification Method:

Protein A purification

CloneNo.:

241402D3

Applications

Tested Applications:

IF/ICC, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using POLR3C antibody (84367-7-RR, Clone: 241402D3) at dilution of 1:400 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 84367-7-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 84367-7-RR against Human POLR3C were performed. The affinity constant is below 1 pM.