

For Research Use Only

Mouse CD28 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84368-6-PBS



Basic Information

Catalog Number: 84368-6-PBS	GenBank Accession Number: NM_007642.4	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 12487	CloneNo.: 241706C11
Source: Rabbit	UNIPROT ID: P31041	
Isotype: IgG	Full Name: CD28 antigen	
	Calculated MW: 25 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
mouse

Product Information

84368-6-PBS targets CD28 as part of a matched antibody pair:

MP01248-4: 84368-7-PBS capture and 84368-6-PBS detection (validated in Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01248-4, Mouse CD28 Recombinant Matched Antibody Pair - PBS only. 84368-7-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg1408. 84368-6-PBS was HRP conjugated as the detection antibody. Range: 78.1-2500 pg/mL.



The mean CD28 concentration was determined to be 105.0 ng/mL in EL-4 cell extract based on a 5.2 mg/mL extract load.