

For Research Use Only

MYH13 Recombinant antibody

Catalog Number: 84379-2-RR



Basic Information

Catalog Number:

84379-2-RR

Size:

100ul , Concentration: 500 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_003802

GeneID (NCBI):

8735

UNIPROT ID:

Q9UIK3

Full Name:

myosin, heavy chain 13, skeletal muscle

Calculated MW:

224 kDa

Observed MW:

230-240 kDa

Purification Method:

Protein A purification

CloneNo.:

241179A3

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse eye tissue, rat eye tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84379-2-RR (MYH13 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.