

For Research Use Only

RNASE3 Recombinant antibody

Catalog Number: 84381-1-RR



Basic Information

Catalog Number:

84381-1-RR

Size:

100ul , Concentration: 1000 µg/ml by 6037

Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_002935

GeneID (NCBI):

6037

UNIPROT ID:

P12724

Full Name:

ribonuclease, RNase A family, 3
(eosinophil cationic protein)

Calculated MW:

18 kDa

Observed MW:

18-20 kDa

Purification Method:

Protein A purification

CloneNo.:

241182E4

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : U-937 cells,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



U-937 cells were subjected to SDS PAGE followed by western blot with 84381-1-RR (RNASE3 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.