

For Research Use Only

Mouse PDGFR beta Recombinant antibody, PBS Only (Detector)

Catalog Number: 84388-1-PBS



Basic Information

Catalog Number:

84388-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_008809.2

GeneID (NCBI):

18596

UNIPROT ID:

P05622-1

Full Name:

platelet derived growth factor receptor, beta polypeptide

Calculated MW:

123kDa

Purification Method:

Protein A purification

CloneNo.:

241754A5

Applications

Tested Applications:

Cytometric bead array, Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

mouse

Product Information

84388-1-PBS targets PDGFR beta as part of a matched antibody pair:

MP01260-1: 84388-3-PBS capture and 84388-1-PBS detection (validated in Cytometric bead array, Sandwich ELISA)

MP01260-2: 84388-4-PBS capture and 84388-1-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Sandwich ELISA standard curve of MP01260-1, Mouse PDGFR beta Recombinant Matched Antibody Pair - PBS only. 84388-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg1533. 84388-1-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL



Serum of six mice was measured. The mouse PDGFR beta concentration of detected samples was determined to be 34.36 ng/mL with a range of 25.40 - 43.60 ng/mL.



Cytometric bead array standard curve of MP01260-1, MOUSE PDGFR beta Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84388-3-PBS. Detection antibody: 84388-1-PBS. Standard: Eg1533. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01260-2, MOUSE PDGFR beta Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84388-4-PBS. Detection antibody: 84388-1-PBS. Standard: Eg1533. Range: 0.625-40 ng/mL