

For Research Use Only

Rat NKp46/NCR1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84416-3-PBS



Basic Information

Catalog Number: 84416-3-PBS	GenBank Accession Number: NM_057199.1	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 117547	CloneNo.: 241811E7
Source: Rabbit	UNIPROT ID: Q9Z0H5	
Isotype: IgG	Full Name: natural cytotoxicity triggering receptor 1	
	Calculated MW: 37kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
rat

Product Information

84416-3-PBS targets NKp46/NCR1 as part of a matched antibody pair:

MP01299-1: 84416-3-PBS capture and 84416-2-PBS detection (validated in Cytometric bead array)

MP01299-2: 84416-3-PBS capture and 84416-4-PBS detection (validated in Cytometric bead array)

MP01299-3: 84416-3-PBS capture and 84416-1-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP01299-1, RAT NKp46/NCR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84416-3-PBS. Detection antibody: 84416-2-PBS. Standard: Eg1469. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01299-2, RAT NKp46/NCR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84416-3-PBS. Detection antibody: 84416-4-PBS. Standard: Eg1469. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01299-3, RAT NKp46/NCR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84416-3-PBS. Detection antibody: 84416-1-PBS. Standard: Eg1469. Range: 0.313-40 ng/mL