

For Research Use Only

C1orf170 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84425-1-PBS



Basic Information

Catalog Number: 84425-1-PBS	GenBank Accession Number: BC006300	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 84808	CloneNo.: 241613A4
Source: Rabbit	UNIPROT ID: Q5SV97	
Isotype: IgG	Full Name: chromosome 1 open reading frame 170	
Immunogen Catalog Number: AG23833		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

84425-1-PBS targets C1orf170 as part of a matched antibody pair:

MP01296-1: 84425-1-PBS capture and 84425-3-PBS detection (validated in Cytometric bead array)

MP01296-2: 84425-1-PBS capture and 84425-2-PBS detection (validated in Cytometric bead array)

MP01296-3: 84425-1-PBS capture and 84425-4-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP01296-1, C1orf170 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84425-1-PBS. Detection antibody: 84425-3-PBS. Standard: Ag23833. Range: 0.469-60 ng/mL



Cytometric bead array standard curve of MP01296-2, C1orf170 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84425-1-PBS. Detection antibody: 84425-2-PBS. Standard: Ag23833. Range: 0.469-60 ng/mL



Cytometric bead array standard curve of MP01296-3, C1orf170 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84425-1-PBS. Detection antibody: 84425-4-PBS. Standard: Ag23833. Range: 0.469-60 ng/mL