

For Research Use Only

Mouse TNFSF11 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84431-2-PBS



Basic Information

Catalog Number:

84431-2-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_011613.3

GeneID (NCBI):

21943

UNIPROT ID:

O35235

Full Name:

tumor necrosis factor (ligand)
superfamily, member 11

Calculated MW:

35KD

Purification Method:

Protein A purification

CloneNo.:

241791A10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Product Information

84431-2-PBS targets TNFSF11 as part of a matched antibody pair:

MP01276-3: 84431-1-PBS capture and 84431-2-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01276-3, MOUSE TNFSF11/RANKL Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84431-1-PBS. Detection antibody: 84431-2-PBS. Standard: Eg1956. Range: 0.313-40 ng/mL