

For Research Use Only

# Nectin-2/CD112 Recombinant antibody

Catalog Number: 84439-2-RR



## Basic Information

Catalog Number:

84439-2-RR

Size:

100ul

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_008990.3

GeneID (NCBI):

19294

UNIPROT ID:

P32507-1

Full Name:

poliovirus receptor-related 2

Calculated MW:

57kDa

Purification Method:

Protein A purification

CloneNo.:

241798A7

## Applications

Tested Applications:

ELISA

Species Specificity:

mouse

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data