

For Research Use Only

OPA1 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 84457-3-PBS



Basic Information

Catalog Number: 84457-3-PBS	GenBank Accession Number: BC075805	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 4976	CloneNo.: 260763B7
Source: Rabbit	UNIPROT ID: O60313	
Isotype: IgG	Full Name: optic atrophy 1 (autosomal dominant)	
Immunogen Catalog Number: AG26887	Calculated MW: 960 aa, 112 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

84457-3-PBS targets OPA1 as part of a matched antibody pair:

MP03472-1: 84457-4-PBS capture and 84457-3-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean OPA1 concentration was determined to be 11.3 ng/mL in A431 cell extract based on a 1.2 mg/mL extract load, 4.2 ng/mL in SKOV-3 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP03472-1, OPA1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84457-4-PBS. Detection antibody: 84457-3-PBS. Standard: Ag26887. Range: 0.156-20 ng/mL.