

For Research Use Only

Mouse IL-3RB/CD131 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84513-2-PBS



Basic Information

Catalog Number:

84513-2-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2173

GenBank Accession Number:

NM_007780.4

GeneID (NCBI):

12983

UNIPROT ID:

P26955

Full Name:

colony stimulating factor 2 receptor, beta, low-affinity (granulocyte-macrophage)

Calculated MW:

99kDa

Purification Method:

Protein A purification

CloneNo.:

241954C11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Product Information

84513-2-PBS targets IL-3RB/CD131 as part of a matched antibody pair:

MP01393-2: 84513-2-PBS capture and 84513-3-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01393-2, MOUSE IL-3RB/CD131 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84513-2-PBS. Detection antibody: 84513-3-PBS. Standard: Eg2173. Range: 0.313-40 ng/mL