

For Research Use Only

Mouse CD64 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84518-3-PBS



Basic Information

Catalog Number: 84518-3-PBS	GenBank Accession Number: NM_010186.5	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 14129	CloneNo.: 241809G4
Source: Rabbit	UNIPROT ID: P26151	
Isotype: IgG	Full Name: Fc receptor, IgG, high affinity I	
	Calculated MW: 45 kDa	
	Observed MW: 68-72 kDa	

Applications

Tested Applications:
WB, IF/ICC, Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Product Information

84518-3-PBS targets CD64 as part of a matched antibody pair:

MP01366-2: 84518-3-PBS capture and 84518-1-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Background Information

Fcγ receptor comprise a multigene family of integral membrane glycoproteins that exhibit complex activation or inhibitory effects on cell functions after aggregation by complexed immunoglobulin G (IgG) (PMID: 17005690). CD64, also known as Fcγ R1A, is a high-affinity receptor for the Fc region of IgG. It is expressed by monocytes/macrophages, activated neutrophils, dendritic cells, and early myeloid cells (PMID: 23293080; 19642859; 7680917). CD64 functions in both innate and adaptive immune responses.


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only



For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84518-3-RR (CD64 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 84518-3-PBS in a different storage buffer formulation.



Cytometric bead array standard curve of MP01366-2, MOUSE CD64 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84518-3-PBS. Detection antibody: 84518-1-PBS. Standard: Eg2107. Range: 0.625-40 ng/mL.