

For Research Use Only

CCR8 Recombinant antibody, PBS Only

Catalog Number: 84545-1-PBS



Basic Information

Catalog Number:

84545-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC074743

GeneID (NCBI):

1237

UNIPROT ID:

P51685

Full Name:

chemokine (C-C motif) receptor 8

Calculated MW:

355 aa, 41 kDa

Observed MW:

41 kDa

Purification Method:

Protein A purification

CloneNo.:

241663E4

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 84545-1-RR (CCR8 antibody) at dilution of 1:800 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 84545-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 84545-1-RR (CCR8 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 84545-1-PBS in a different storage buffer formulation.