

For Research Use Only

IL15RA Recombinant antibody, PBS Only (Detector)

Catalog Number: 84582-3-PBS



Basic Information

Catalog Number: 84582-3-PBS	GenBank Accession Number: BC107777	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 3601	CloneNo.: 241756B10
Source: Rabbit	UNIPROT ID: Q13261	
Isotype: IgG	Full Name: interleukin 15 receptor, alpha	
	Calculated MW: 231 aa, 24 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

84582-3-PBS targets IL15RA as part of a matched antibody pair:

MP01405-2: 84582-1-PBS capture and 84582-3-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01405-2, IL-15RA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84582-1-PBS. Detection antibody: 84582-3-PBS. Standard: RP00015. Range: 0.469-60 ng/mL