

For Research Use Only

# EDDM3B Recombinant antibody, PBS Only (Detector)

Catalog Number: 84597-1-PBS



## Basic Information

**Catalog Number:**

84597-1-PBS

**Size:**

100ug, Concentration: 1 mg/ml by Nanodrop;

**Source:**

Rabbit

**Isotype:**

IgG

**GenBank Accession Number:**

NM\_022360.5

**GeneID (NCBI):**

64184

**UNIPROT ID:**

P56851

**Full Name:**

family with sequence similarity 12, member B (epididymal)

**Calculated MW:**

18kDa

**Purification Method:**

Protein A purification

**CloneNo.:**

241998B5

## Applications

**Tested Applications:**

Cytometric bead array, Indirect ELISA

**Species Specificity:**

human

## Product Information

84597-1-PBS targets EDDM3B as part of a matched antibody pair:

MP01431-1: 84597-2-PBS capture and 84597-1-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage

**Storage:**

Store at -80°C.

**Storage Buffer:**

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01431-1, EDDM3B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84597-2-PBS. Detection antibody: 84597-1-PBS. Standard: Eg2443. Range: 0.469-60 ng/mL