

For Research Use Only

C11orf68 Recombinant antibody

Catalog Number: 84725-4-RR



Basic Information

Catalog Number: 84725-4-RR	GenBank Accession Number: BC010512	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 83638	CloneNo.: 242221F8
Source: Rabbit	UNIPROT ID: Q9H3H3	Recommended Dilutions: WB 1:2000-1:12000
Isotype: IgG	Full Name: chromosome 11 open reading frame 68	
Immunogen Catalog Number: AG36023	Observed MW: 35 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : A431 cells, mouse liver tissue, rat pancreas tissue
Species Specificity: human, mouse, rat	

Background Information

BLES03 (basophilic leukemia-expressed protein 03) is a human protein of unknown function. Recent studies have shown that BLES03 is a novel eukaryotic phosphothreonine cleavage enzyme involved in the biosynthesis of dehydrogenated amino acids. The crystal structure of BLES03 has been resolved by single-wavelength anomalous scattering (SAD) with a resolution of 2.5 Å.

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84725-4-RR (C11orf68 antibody) at dilution of 1:6000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 84725-4-RR against Human C11orf68 were performed. The affinity constant is 22.8 nM.