

For Research Use Only

Mono-Methyl-Histone H3 (Lys4) Recombinant antibody, PBS Only

Catalog Number: 84740-1-PBS



Basic Information

Catalog Number:

84740-1-PBS

Size:

100ug , Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC066245

GeneID (NCBI):

8350

UNIPROT ID:

P68431

Full Name:

histone cluster 1, H3a

Observed MW:

15 kDa

Purification Method:

Protein A purification

CloneNo.:

241940G8

Applications

Tested Applications:

WB, Dot Blot, Indirect ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Dot blot analysis was used to confirm the specificity of 84740-1-RR Mono-Methyl-Histone H3 (Lys4) antibody. Peptides were spotted onto NC and probed with antibody at 1 µg/mL. The amount of peptide (ug/mL) spotted is indicated next to each row. This data was developed using the same antibody clone with 84740-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 84740-1-RR (Mono-Methyl-Histone H3 (Lys4) antibody) at dilution of 1:15000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 84740-1-PBS in a different storage buffer formulation.