

For Research Use Only

SLC7A7 Recombinant antibody, PBS Only

Catalog Number: 84743-1-PBS



Basic Information

Catalog Number:

84743-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC010107

GeneID (NCBI):

9056

UNIPROT ID:

Q9UM01

Full Name:

solute carrier family 7 (cationic amino acid transporter, y+ system), member 7

Calculated MW:

511 aa, 56 kDa

Observed MW:

50-60 kDa

Purification Method:

Protein A purification

CloneNo.:

241946F7

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84743-1-RR (SLC7A7 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 84743-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 84743-1-RR (SLC7A7 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 84743-1-PBS in a different storage buffer formulation.