

For Research Use Only

SPAG4-Specific Recombinant antibody, PBS Only

Catalog Number: 84842-1-PBS



Basic Information

Catalog Number:

84842-1-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_003116

GeneID (NCBI):

6676

UNIPROT ID:

Q9NPE6

Full Name:

sperm associated antigen 4

Calculated MW:

48 kDa

Purification Method:

Protein A purification

CloneNo.:

241920E1

Applications

Tested Applications:

IHC, Indirect ELISA

Species Specificity:

human, mouse

Background Information

SPAG4 may assist the organization and assembly of outer dense fibers (ODFs), a specific structure of the sperm tail.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded mouse testis tissue slide using 84842-1-RR (SPAG4-Specific antibody) at dilution of 1:400 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 84842-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded mouse testis tissue slide using 84842-1-RR (SPAG4-Specific antibody) at dilution of 1:400 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 84842-1-PBS in a different storage buffer formulation.