

For Research Use Only

# IL-10RB Recombinant antibody, PBS Only (Capture)

Catalog Number: 84858-2-PBS



## Basic Information

<b>Catalog Number:</b> 84858-2-PBS	<b>GenBank Accession Number:</b> BC001903	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ug , Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 3588	<b>CloneNo.:</b> 242302B1
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q08334	
<b>Isotype:</b> IgG	<b>Full Name:</b> interleukin 10 receptor, beta	
	<b>Calculated MW:</b> 37 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Sandwich ELISA, Indirect ELISA, Sample test

**Species Specificity:**  
human

## Product Information

84858-2-PBS targets IL-10RB as part of a matched antibody pair:

MP01631-1: 84858-2-PBS capture and 84858-1-PBS detection (validated in Cytometric bead array, Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.


## Storage

**Storage:**  
Store at -80°C.  
**Storage Buffer:**  
PBS Only


For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


# Selected Validation Data



Sandwich ELISA standard curve of MP01631-1, Human IL-10RB Recombinant Matched Antibody Pair - PBS only. 84858-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg3446. 84858-1-PBS was HRP conjugated as the detection antibody. Range: 0.313-20 ng/mL



The mean IL-10RB concentration was determined to be 5.84 ng/mL in THP-1 cell extract based on a 4.40 mg/mL extract load.



Cytometric bead array standard curve of MP01631-1, IL-10RB Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84858-2-PBS. Detection antibody: 84858-1-PBS. Standard: Eg3446. Range: 0.156-20 ng/mL