

For Research Use Only

PMS2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84894-7-PBS



Basic Information

Catalog Number: 84894-7-PBS	GenBank Accession Number: BC093921	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 5395	CloneNo.: 242455C9
Source: Rabbit	UNIPROT ID: P54278	
Isotype: IgG	Full Name: PMS2 postmeiotic segregation increased 2 (S. cerevisiae)	
Immunogen Catalog Number: AG33508	Calculated MW: 862 aa, 96 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Product Information

84894-7-PBS targets PMS2 as part of a matched antibody pair:

MP01663-3: 84894-6-PBS capture and 84894-7-PBS detection (validated in Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01663-3, Human PMS2 Recombinant Matched Antibody Pair - PBS only. 84894-6-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag33508. 84894-7-PBS was HRP conjugated as the detection antibody. Range: 0.391-25 ng/mL.



The mean PMS2 concentration was determined to be 32.56 ng/mL in HT-29 cell extract based on a 1.50 mg/mL extract load, 67.06 ng/mL in HCT 116 cell extract based on a 1.40 mg/mL extract load and 64.46 ng/mL in U-251 cell extract based on a 1.40 mg/mL extract load.