

For Research Use Only

# TRIM56 Recombinant antibody, PBS Only

Catalog Number: 84896-7-PBS



## Basic Information

<b>Catalog Number:</b> 84896-7-PBS	<b>GenBank Accession Number:</b> BC048194	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ug , Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 81844	<b>CloneNo.:</b> 242497D3
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q9BRZ2	
<b>Isotype:</b> IgG	<b>Full Name:</b> tripartite motif-containing 56	
<b>Immunogen Catalog Number:</b> AG21992	<b>Calculated MW:</b> 755 aa, 81 kDa	
	<b>Observed MW:</b> 81 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA

**Species Specificity:**  
human

## Storage


**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed MCF-7 cells using TRIM56 antibody (84896-7-RR, Clone: 242497D3 ) at dilution of 1:200 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 84896-7-PBS in a different storage buffer formulation.