

For Research Use Only

TRIM56 Recombinant antibody

Catalog Number: 84896-7-RR



Basic Information

Catalog Number:

84896-7-RR

Size:

100ul , Concentration: 1000 µg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21992

GenBank Accession Number:

BC048194

GeneID (NCBI):

81844

UNIPROT ID:

Q9BRZ2

Full Name:

tripartite motif-containing 56

Calculated MW:

755 aa, 81 kDa

Observed MW:

81 kDa

Purification Method:

Protein A purification

CloneNo.:

242497D3

Applications

Tested Applications:

ELISA

Species Specificity:

human

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed MCF-7 cells using TRIM56 antibody (84896-7-RR, Clone: 242497D3) at dilution of 1:200 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2).