

For Research Use Only

PSPH Recombinant antibody, PBS Only

Catalog Number: 84899-5-PBS



Basic Information

Catalog Number:

84899-5-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5972

GenBank Accession Number:

BC063614

GeneID (NCBI):

5723

UNIPROT ID:

P78330

Full Name:

phosphoserine phosphatase

Calculated MW:

25 kDa

Observed MW:

25-28 kDa

Purification Method:

Protein A purification

CloneNo.:

242518C4

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84899-5-RR (PSPH antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 84899-5-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 84899-5-RR against Human PSPH were performed. The affinity constant is below 1 pM.