

For Research Use Only

TCOF1 Recombinant antibody, PBS Only

Catalog Number: 84914-5-PBS



Basic Information

Catalog Number: 84914-5-PBS	GenBank Accession Number: BC014559	Purification Method: Protein A purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 6949	CloneNo.: 242201F12
Source: Rabbit	UNIPROT ID: Q13428	
Isotype: IgG	Full Name: Treacher Collins-Franceschetti syndrome 1	
Immunogen Catalog Number: AG1447	Calculated MW: 97 kDa, 144 kDa, 152 kDa	

Applications

Tested Applications:
IF/ICC, Indirect ELISA

Species Specificity:
human

Background Information

TCOF1 also known as treacle, is a nucleolar protein with a LIS1 homology domain. TCOF1 is involved in ribosomal DNA gene transcription by interacting with upstream binding factor (UBF). Mutations in TCOF1 have been associated with Treacher Collins syndrome, a disorder which includes abnormal craniofacial development(PMID: 31218872)

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed MCF-7 cells using TCOF1 antibody (84914-5-RR, Clone: 242201F12) at dilution of 1:500 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 84914-5-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 84914-5-RR against Human TCOF1 were performed. The affinity constant is 8.67 nM.