

For Research Use Only

TRIM27 Recombinant antibody

Catalog Number: 84927-5-RR



Basic Information

Catalog Number:

84927-5-RR

Size:

100ul, Concentration: 1000 µg/ml by

Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2849

GenBank Accession Number:

BC013580

GeneID (NCBI):

5987

UNIPROT ID:

P14373

Full Name:

tripartite motif-containing 27

Calculated MW:

513 aa, 58 kDa

Observed MW:

58 kDa

Purification Method:

Protein A purification

CloneNo.:

242538B5

Recommended Dilutions:

WB 1:2000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HEK-293 cells, MCF-7 cells, HeLa cells, Raji cells, Jurkat cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84927-5-RR (TRIM27 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.