

For Research Use Only

# GOLGA3 Recombinant antibody

Catalog Number: 84955-4-RR



## Basic Information

Catalog Number:

84955-4-RR

Size:

100ul , Concentration: 1000 µg/ml by 2802

Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15531

GenBank Accession Number:

BC142658

GeneID (NCBI):

2802

UNIPROT ID:

Q08378

Full Name:

golgi autoantigen, golgin subfamily a, 3

Calculated MW:

1498 aa, 167 kDa

Observed MW:

167 kDa

Purification Method:

Protein A purification

CloneNo.:

242603C7

Recommended Dilutions:

WB 1:1000-1:8000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HEK-293 cells, HeLa cells, L02 cells, PC-3 cells

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84955-4-RR (GOLGA3 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 84955-4-RR against Human GOLGA3 were performed. The affinity constant is 0.475 nM.