

For Research Use Only

TFB1M Recombinant antibody

Catalog Number: 84957-4-RR



Basic Information

Catalog Number:

84957-4-RR

Size:

100ul , Concentration: 1000 µg/ml by

Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9785

GenBank Accession Number:

BC017788

GeneID (NCBI):

51106

UNIPROT ID:

Q8WVM0

Full Name:

transcription factor B1, mitochondrial

Calculated MW:

346 aa, 40 kDa

Observed MW:

40 kDa

Purification Method:

Protein A purification

CloneNo.:

242524A10

Recommended Dilutions:

WB 1:2000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HEK-293T cells, HepG2 cells, K-562 cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84957-4-RR (TFB1M antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 84957-4-RR against Human TFB1M were performed. The affinity constant is 2.99 nM.