

For Research Use Only

SUGT1 Recombinant antibody

Catalog Number: 84963-5-RR



Basic Information

Catalog Number:

84963-5-RR

Size:

100ul , Concentration: 1000 µg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2275

GenBank Accession Number:

BC000911

GeneID (NCBI):

10910

UNIPROT ID:

Q9Y2Z0

Full Name:

SGT1, suppressor of G2 allele of SKP1 (S. cerevisiae)

Calculated MW:

41 kDa

Observed MW:

38-41 kDa

Purification Method:

Protein A purification

CloneNo.:

242590G3

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HCT 116 cells, HeLa cells, HEK-293 cells, Jurkat cells, K-562 cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84963-5-RR (SUGT1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.