

For Research Use Only

MAS1L Recombinant antibody

Catalog Number: 85003-1-RR



Basic Information

Catalog Number:

85003-1-RR

Size:

100ul , Concentration: 1000 µg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_052967

GeneID (NCBI):

116511

UNIPROT ID:

P35410

Full Name:

MAS1 oncogene-like

Calculated MW:

42 kDa

Observed MW:

42 kDa

Purification Method:

Protein A purification

CloneNo.:

242360D1

Recommended Dilutions:

WB 1:2000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : HeLa cells, HepG2 cells, HEK-293 cells, COLO 320 cells, NIH/3T3 cells

Background Information

MAS1L, also named MRG, belongs to the G-protein coupled receptor 1 family and Mas subfamily.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85003-1-RR (MAS1L antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.