

For Research Use Only

# Mouse LAIR1 Recombinant antibody, PBS Only

Catalog Number: 85010-2-PBS



## Basic Information

Catalog Number:

85010-2-PBS

Size:

100ug, Concentration: 1 mg/ml by  
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG1879

GenBank Accession Number:

NM\_001113474.1

GeneID (NCBI):

52855

UNIPROT ID:

Q8BG84-1

Full Name:

leukocyte-associated Ig-like receptor

1

Calculated MW:

30kDa

Purification Method:

Protein A purification

CloneNo.:

242226E5

## Applications

Tested Applications:

Indirect ELISA

Species Specificity:

mouse

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech  
Group brand and is not available to purchase from any  
other manufacturer.

