

For Research Use Only

PDLIM1,CLP36 Recombinant antibody, PBS Only (Detector)

Catalog Number:85027-1-PBS



Basic Information

Catalog Number: 85027-1-PBS	GenBank Accession Number: BC000915	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 9124	CloneNo.: 242585C8
Source: Rabbit	UNIPROT ID: O00151	
Isotype: IgG	Full Name: PDZ and LIM domain 1	
Immunogen Catalog Number: AG2274	Calculated MW: 329 aa, 36 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

85027-1-PBS targets PDLIM1,CLP36 as part of a matched antibody pair:

MP01749-1: 85027-3-PBS capture and 85027-1-PBS detection (validated in Cytometric bead array)

MP01749-2: 85027-2-PBS capture and 85027-1-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01749-1, PDLIM1/CLP36 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85027-3-PBS. Detection antibody: 85027-1-PBS. Standard: Ag2274. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01749-2, PDLIM1/CLP36 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85027-2-PBS. Detection antibody: 85027-1-PBS. Standard: Ag2274. Range: 0.313-40 ng/mL