

For Research Use Only

LAMB3 Recombinant antibody

Catalog Number: 85055-4-RR



Basic Information

Catalog Number:

85055-4-RR

Size:

100ul , Concentration: 1000 ug/ml by

Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25192

GenBank Accession Number:

BC075838

GeneID (NCBI):

3914

UNIPROT ID:

Q13751

Full Name:

laminin, beta 3

Calculated MW:

1172 aa, 130 kDa

Observed MW:

130-140 kDa

Purification Method:

Protein A purification

CloneNo.:

242442A8

Recommended Dilutions:

WB 1:2000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, rat

Positive Controls:

WB : A431 cells, U2OS cells, rat lung tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85055-4-RR (LAMB3 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 85055-4-RR against Human LAMB3 were performed. The affinity constant is 0.244 nM.