

For Research Use Only

ADSS Recombinant antibody

Catalog Number: 85056-3-RR



Basic Information

Catalog Number:

85056-3-RR

Size:

100ul, Concentration: 1000 µg/ml by

Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9298

GenBank Accession Number:

BC012356

GeneID (NCBI):

159

UNIPROT ID:

P30520

Full Name:

adenylosuccinate synthase

Calculated MW:

456 aa, 50 kDa

Observed MW:

45 kDa

Purification Method:

Protein A purification

CloneNo.:

242492F12

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HeLa cells, HepG2 cells, HEK-293T cells, K-562 cells, C6 cells, NIH/3T3 cells, mouse skeletal muscle tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85056-3-RR (ADSS antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.