

For Research Use Only

THRAP3 Recombinant antibody, PBS Only

Catalog Number: 85093-1-PBS



Basic Information

Catalog Number:

85093-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_005119

GeneID (NCBI):

9967

UNIPROT ID:

Q9Y2W1

Full Name:

thyroid hormone receptor associated protein 3

Calculated MW:

109 kDa

Observed MW:

109-130 kDa

Purification Method:

Protein A purification

CloneNo.:

242268B5

Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human

Background Information

THRAP3, also named as TRAP150, plays a role in transcriptional coactivation.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85093-1-RR (THRAP3 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85093-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using THRAP3 antibody (85093-1-RR, Clone: 242268B5) at dilution of 1:600 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 85093-1-PBS in a different storage buffer formulation.