

For Research Use Only

ACSS1 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 85102-3-PBS



Basic Information

Catalog Number: 85102-3-PBS	GenBank Accession Number: BC039261	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 84532	CloneNo.: 242361F8
Source: Rabbit	UNIPROT ID: Q9NUB1	
Isotype: IgG	Full Name: acyl-CoA synthetase short-chain family member 1	
Immunogen Catalog Number: AG10853	Calculated MW: 689 aa, 75 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

85102-3-PBS targets ACSS1 as part of a matched antibody pair:

MP01848-1: 85102-2-PBS capture and 85102-3-PBS detection (validated in Cytometric bead array)

MP01848-2: 85102-3-PBS capture and 85102-1-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01848-2, ACS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85102-3-PBS. Detection antibody: 85102-1-PBS. Standard: Ag10853. Range: 3.125-200 ng/mL.



Cytometric bead array standard curve of MP01848-1, ACS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85102-2-PBS. Detection antibody: 85102-3-PBS. Standard: Ag10853. Range: 1.563-100 ng/mL.