

For Research Use Only

Mouse Foxp3 Recombinant antibody, PBS Only

Catalog Number: 85128-4-PBS



Basic Information

Catalog Number:

85128-4-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

GeneID (NCBI):

20371

Full Name:

forkhead box P3

Observed MW:

45 kDa

Purification Method:

Protein A purification

CloneNo.:

242674B10

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Background Information

FOXP3, also named as IPEX, JM2 and Scurfin, is a transcription factor. FOXP3 plays a critical role in the control of immune response. Defects in FOXP3 are the cause of immunodeficiency polyendocrinopathy, enteropathy, X-linked syndrome (IPEX) which also known as X-linked autoimmunity-immunodeficiency syndrome.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85128-4-RR (Foxp3 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85128-4-PBS in a different storage buffer formulation.