

For Research Use Only

CARM1 Recombinant antibody

Catalog Number: 85146-3-RR



Basic Information

Catalog Number:

85146-3-RR

Size:

100ul , Concentration: 1000 µg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_199141

GeneID (NCBI):

10498

UNIPROT ID:

Q86X55

Full Name:

coactivator-associated arginine
methyltransferase 1

Calculated MW:

66 kDa

Observed MW:

63-65 kDa

Purification Method:

Protein A purification

CloneNo.:

242651F10

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HeLa cells, HEK-293T cells, HCT 116 cells, MCF-7
cells, NCI-H1299 cells, NIH/3T3 cells, PC-12 cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85146-3-RR (CARM1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.