

For Research Use Only

MAGEB3 Recombinant antibody, PBS Only

Catalog Number: 85191-1-PBS



Basic Information

Catalog Number: 85191-1-PBS	GenBank Accession Number: BC074756	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 4114	CloneNo.: 242499A11
Source: Rabbit	UNIPROT ID: O15480	
Isotype: IgG	Full Name: melanoma antigen family B, 3	
Immunogen Catalog Number: AG12476	Calculated MW: 346 aa, 39 kDa	
	Observed MW: 39 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, mouse

Background Information

MAGE genes encode for tumor-rejection antigens and are expressed in tumors of different histologic types and Testis/placenta, but not in other normal tissues. MAGEB3, also named CT3.5, is a member of the MAGEB gene family.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse testis tissue was subjected to SDS PAGE followed by western blot with 85191-1-RR (MAGEB3 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85191-1-PBS in a different storage buffer formulation.