

For Research Use Only

OAT Recombinant antibody

Catalog Number: 85218-2-RR



Basic Information

Catalog Number:

85218-2-RR

Size:

100ul , Concentration: 1000 µg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9981

GenBank Accession Number:

BC000964

GeneID (NCBI):

4942

UNIPROT ID:

P04181

Full Name:

ornithine aminotransferase (gyrate atrophy)

Calculated MW:

48 kDa

Observed MW:

48 kDa

Purification Method:

Protein A purification

CloneNo.:

242843D8

Recommended Dilutions:

WB 1:1000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : K-562 cells, HCT 116 cells, HepG2 cells, mouse colon tissue, mouse lung tissue, rat colon tissue, rat lung tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85218-2-RR (OAT antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.