

For Research Use Only

Phospho-IRS1 (Ser307) Recombinant antibody, PBS Only

Catalog Number: 85238-1-PBS



Basic Information

Catalog Number:

85238-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC053895

GeneID (NCBI):

3667

UNIPROT ID:

P35568

Full Name:

insulin receptor substrate 1

Calculated MW:

1242 aa, 132 kDa

Observed MW:

180 kDa

Purification Method:

Protein A purification

CloneNo.:

242124E4

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Background Information

Serine phosphorylation of insulin receptor substrate-1 (IRS-1) inhibits insulin signal transduction in a variety of cell backgrounds, which might contribute to peripheral insulin resistance. However, because of the large number of potential phosphorylation sites, the mechanism of inhibition has been difficult to determine. It has been confirmed that JNK and IKK are involved in the phosphorylation of IRS-1 (Ser307/312). Additionally, activation of Akt, ERK, mTOR, and PKC ζ correlates to IRS-1 phosphorylation at Ser307/312, Ser267/270, and Ser612/616 in rodent/human IRS-1. (PMID: 11606564, PMID: 12714600)

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Non-treated NIH/3T3 cells and TNF alpha treated NIH/3T3 cells were subjected to SDS PAGE followed by western blot with 85238-1-RR (Phospho-IRS1 (Ser307) antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. The membrane was stripped and re-blotted with Lamin B1 (66095-1-Ig) antibody as a loading control. This data was developed using the same antibody clone with 85238-1-PBS in a different storage buffer formulation.



Non-treated and TNF alpha treated HeLa cells were subjected to SDS PAGE followed by western blot with 85238-1-RR (Phospho-IRS1 (Ser307) antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. The membrane was stripped and re-blotted with Vinculin (66305-1-Ig) antibody as a loading control. This data was developed using the same antibody clone with 85238-1-PBS in a different storage buffer formulation.