

For Research Use Only

# IRF1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85254-5-PBS



## Basic Information

<b>Catalog Number:</b> 85254-5-PBS	<b>GenBank Accession Number:</b> BC009483	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ug , Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 3659	<b>CloneNo.:</b> 250167G11
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P10914	
<b>Isotype:</b> IgG	<b>Full Name:</b> IRF 1	
<b>Immunogen Catalog Number:</b> AG1879	<b>Calculated MW:</b> 325 aa, 37 kDa	

## Applications

**Tested Applications:**  
Sandwich ELISA, Indirect ELISA, Sample test

**Species Specificity:**  
human

## Product Information

85254-5-PBS targets IRF1 as part of a matched antibody pair.

MP02217-1: 85254-5-PBS capture and 85254-4-PBS detection (validated in Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.


## Storage

**Storage:**  
Store at -80°C.  
**Storage Buffer:**  
PBS only, pH7.3


For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP02217-1, Human IRF1 Recombinant Matched Antibody Pair - PBS only. 85254-5-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag1879. 85254-4-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL.



The mean IRF1 concentration was determined to be 12.24 ng/mL in HL-60 cell extract based on a 1.50 mg/mL extract load and 8.45 ng/mL in THP-1 cell extract based on a 3.40 mg/mL extract load.