

For Research Use Only

ART4 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85273-5-PBS



Basic Information

Catalog Number: 85273-5-PBS	GenBank Accession Number: NM_021071.4	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 420	CloneNo.: 250103C9
Source: Rabbit	UNIPROT ID: Q93070	
Isotype: IgG	Full Name: ADP-ribosyltransferase 4 (Dombrock blood group)	
Immunogen Catalog Number: EG2665	Calculated MW: 36 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Product Information

85273-5-PBS targets ART4 as part of a matched antibody pair:

MP02589-1: 85273-5-PBS capture and 85273-4-PBS detection (validated in Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02589-1, Human ART4/CD297 Recombinant Matched Antibody Pair - PBS only. 85273-5-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg2665. 85273-4-PBS was HRP conjugated as the detection antibody. Range: 0.781-50 ng/mL



The mean ART4/CD297 concentration was determined to be 251.26 ng/mL in HeLa cell extract based on a 1.20 mg/mL extract load and 280.28 ng/mL in Jurkat cell extract based on a 1.20 mg/mL extract load.